

PETER H. CAIN, D.D.S.
Orthodontics and Dentofacial Orthopedics
5 DURHAM ROAD
GUILFORD, CONNECTICUT 06437
—
Telephone (203) 453-5955
Facsimile (203) 458-0426
shorelineorthoct.com

RED, WHITE & BLUE INVISIBLE APPLIANCES

Over the course of your treatment you will be receiving several invisible appliances which have been specially produced from models of your teeth – they only fit you. The purpose, or function, of the appliances is to achieve alignment, or straightening, of your teeth.

The appliances are activated each time you firmly seat them over your teeth. The special material has been pre-programmed to coax you teeth toward a more ideal and pleasing position. While the amount of movement each appliance can produce is measured in fractions of a millimeter, the force will be noticeable and may result in some temporary minor discomfort. After approximately 4-6 weeks your teeth should have moved to the designed correction. This is your indication to move on to the next unit.

Your appliances should be worn at all times – even when you sleep. You should remove the appliance when eating. If you lose an appliance, you need to call us immediately. Lost wearing time can extend treatment time and compromise results.

Clean your retainer prior to each insertion with a soft bristle toothbrush using water or a small amount of toothpaste. Rinse thoroughly with water after cleaning. Do not use denture cleaners to clean your retainer or soak it in mouthwash. These products can damage the surface of the retainer, causing it to become dull and more visible.

Keep your retainer away from hot water, any heat source, or harsh chemicals.

Keep your retainer out of reach of young children and pets.

Remove your retainer for eating and drinking, except when drinking water.

REMEMBER, if the appliance isn't in your mouth it should be in the carrying case to prevent its being accidentally lost.

Great results are in your hands and mouth. If your situation requires additional appliances, it is important to wear your final appliance until we supply you with a new series.